Pressure gauge for sanitary applications
Models PG43SA-S, PG43SA-C

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© 12/2011 WIKA Alexander Wiegand SE & Co. KG

All rights reserved.

WIKA® is a registered trademark in various countries.

Prior to starting any work, read the operating instructions!

Keep for later use!

Contents

1. General information	4
2. Design and function	6
3. Safety	9
4. Transport, packaging and storage	11
5. Commissioning, operation	12
6. Faults	17
7. Maintenance and cleaning	17
8. Dismounting, return and disposal	19
9. Specifications	20

Declarations of conformity can be found online at www.wika.com.

1. General information

Supplementary documentation:

- ▶ Please follow all the documentation included within the scope of delivery.

With versions for hazardous areas, also follow the additional operating instructions

1. General information

- The instrument described in the operating instructions has been designed and manufactured using state-of-the-art technology. All components are subject to stringent quality and environmental criteria during production. Our management systems are certified in accordance with ISO 9001 and ISO 14001.
- These operating instructions contain important information on handling the instrument. Working safely requires that all safety notes and work instructions are observed.
- Unauthorised modifications to the product will invalidate the warranty and result in loss of services.
- Observe the relevant local accident prevention regulations and general safety regulations for the instrument's range of use.
- The operating instructions are part of the product and must be kept in the immediate vicinity of the instrument and readily accessible to skilled personnel at any time. Pass the operating instructions on to the next operator or owner of the instrument.
- Skilled personnel must have carefully read and understood the operating instructions prior to beginning any work.
- If available, the provided supplier documentation is also considered to be part of the product in addition to these operating instructions.
- The general terms and conditions contained in the sales documentation shall apply.
- Subject to technical modifications.

2. Design and function

EN

Notes per Pressure Equipment Directive

The instruments are defined as “pressure accessories” per article 1, para. 2.5.

The measuring characteristics are in accordance with the EN 837-3 standard.

■ Further information:

- Internet address: www.wika.de / www.wika.com
- Relevant data sheets: PM 04.15, PM 04.16
- Application consultant: Tel.: +49 9372 132-0
info@wika.de

1.1 Abbreviations, definitions

■ Bulleted list

► Instruction

1. ... x. Follow the instruction step by step

1.2 Explanation of symbols

WARNING!

... indicates a potentially dangerous situation that can result in serious injury or death, if not avoided.

NOTE!

... indicates a potentially dangerous situation that can result in damage to property or the environment, if not avoided.

INFORMATION!

... points out useful tips, recommendations and information for efficient and trouble-free operation.

2. Design and function

EN

2.1 Description

The model PG43SA-S and PG43SA-C pressure gauges are designed especially for the requirements of the food and beverage, pharmaceutical and biotechnology industries.

The pressure gauge is particularly suitable for the special conditions of CIP/SIP cleaning processes, such as chemical stability towards cleaning liquids and high temperatures. The flush metal diaphragm element is directly welded to the process connection. This guarantees a crevice-free connection of process connection and diaphragm element; additional internal sealings are not required.

In special versions, the models PG43SA-S and PG43SA-C are autoclavable.

2.2 Model PG43SA-S with vent hole

In order to avoid any change in pressure within the case, the measuring instrument features a vent hole at the case bottom. In case of condensation formation within the case or ingress of low water quantities from outside, this can be drained through the vent hole at the rear side.


2. Design and function

2.3 Model PG43SA-S with vent valve


Vent valve

Vent valve closed

Vent valve open

EN

For versions with vent valve, the case will be vented to or sealed from the atmosphere by opening or closing the vent valve.

- To close: Turn approx. 3 turns in the direction of the “CLOSE” arrow, until the USIT sealing is pressed together fully.
- To open: Turn approx. 3 turns in the direction of the “OPEN” arrow, until the retaining ring prevents the vent valve from being unscrewed further and thus lost.


For any external cleaning (“wash-down”), e.g. in a cabinet washer, the vent valve must be closed.

2. Design and function

2.4 Model PG43SA-C with vent valve

EN

Only ever operate the vent valve via the spanner flats and to the prescribed torque using an appropriate tool.

By opening the vent valve, the pressure difference between the atmosphere and the sealed case will be eliminated.

Maintaining the ingress protection is only ensured with a closed vent valve.

Only through a brief venting of the case before each pressure measurement can the pressure be read correctly.

2.5 Scope of delivery

Cross-check scope of delivery with delivery note.

3. Safety

3. Safety

EN

3.1 Intended use

The models PG43SA-S and PG43SA-C pressure gauges are used for measuring pressure particularly in sanitary applications.

The instrument has been designed and built solely for the intended use described here, and may only be used accordingly.

The manufacturer shall not be liable for claims of any type based on operation contrary to the intended use.

3.2 Improper use

WARNING!

Injuries through improper use

Improper use of the instrument can lead to hazardous situations and injuries.

- Refrain from unauthorised modifications to the instrument.

Any use beyond or different to the intended use is considered as improper use.

3.3 Compliance with the 3-A Sanitary Standard

For a connection per 3-A Sanitary Standard, the following seals have to be used:

- For milk thread fittings per DIN 11851, suitable profile seals have to be used (e.g. SKS Komponenten BV or Kieselmann GmbH).
- For fittings per IDF, seals with support ring per ISO 2853 have to be used.

Note: Connections per Svensk Standard, APV RJT and NEUMO Connect S are not compliant with the 3-A Sanitary Standard.

3.4 Compliance with EHEDG conformity

To ensure EHEDG-compliant use of the instruments, the general specifications from the EHEDG position paper “Easy Cleanable Pipe Couplings and Process Connections” must be observed.

For an EHEDG compliant connection, seals in accordance with the current EHEDG policy document must be used.

Seals for connections per ISO 2852, DIN 32676 and BS 4825 part 3 are, e.g., manufactured by Combifit International B.V.

A manufacturer of seals for connections per DIN 11851 is, e.g., Kieselmann GmbH. A manufacturer of VARINLINE® seals is, e.g., GEA Tuchenhausen GmbH. A manufacturer of NEUMO BioConnect® seals is, e.g., Neumo GmbH & Co. KG.

3.5 Personnel qualification

WARNING!

Risk of injury should qualification be insufficient

Improper handling can result in considerable injury and damage to property.

- The activities described in these operating instructions may only be carried out by skilled personnel who have the qualifications described below.

Skilled personnel

Skilled personnel, authorised by the operator, are understood to be personnel who, based on their technical training, knowledge of measurement and control technology and on their experience and knowledge of country-specific regulations, current standards and directives, are capable of carrying out the work described and independently recognising potential hazards.

4. Transport, packaging and storage

4.1 Transport

Check the measuring instrument for any damage that may have been caused by transport.

Obvious damage must be reported immediately.

Note!

Damage due to improper transport

With improper transport, a high level of damage to property can occur.

- ▶ When unloading packed goods upon delivery as well as during internal transport, proceed carefully and observe the symbols on the packaging.
- ▶ With internal transport, observe the instructions in chapter 4.2 "Packaging and storage".

4.2 Packaging and storage

Do not remove packaging until just before mounting.

Keep the packaging, especially the diaphragm protection (if available). This will provide optimum protection during transport (e.g. cleaning, change in place of use, sending for repair).

Note!

Damage to the diaphragm during unpacking

With improper unpacking, a high level of damage to property can occur.

- ▶ Avoid touching the diaphragm after removing the diaphragm protection. Impermissible mechanical loading (e.g. pressing the diaphragm) can result in measuring deviations or even in the complete failure of the measuring instrument.

5. Commissioning, operation

EN

The PG43SA-C is optionally packed under sterile conditions, and on the perforation on the rear of the sterile packaging there is a brown-coloured process indicator. By checking this brown colouration you can find out whether the gauge has been packed under sterile conditions.

Store instruments that are packed under sterile conditions in the instrument package (no condensation allowed, storing possible for 6 months).

Permissible conditions at the place of storage

Storage temperature: -40 ... +70 °C [-40 ... +158 °F]

Protect the instrument from dust.

5. Commissioning, operation

5.1 Mounting preparation

- To protect the instrument from mechanical damage keep it in the outer carton until installation.
- The measuring instrument should only be removed from the sterile packaging shortly before mounting at the measuring point.
- With CIP (cleaning in place) we recommend to clean the measuring instrument again after mounting, using the process established at the operator's site.
- With COP (cleaning out of place) we recommend to clean the measuring instrument again before mounting, using the process established at the operator's site.
- On removal from the packaging, and during mounting, particular care must be taken to avoid damage and any mechanical deformation to the diaphragm.
- Model PG43SA-S with vent hole: Check the zero point. In case of a zero offset, carry out a zero point setting. See chapter 5.7 "External zero point setting".

- Model PG43SA-S with vent valve: Open the vent valve for pressure compensation and leave it open during the measurement in order to ensure a correct pressure display. See chapter 2.3 “Model PG43SA-S with vent valve”.
- Model PG43SA-C with vent valve: For pressure compensation, open the vent valve for a short time and then close it again directly. The vent valve remains closed during operation. See chapter 2.4 “Model PG43SA-C with vent valve”.

5.2 Installation

- Do not damage the diaphragm; scratches on the diaphragm (e.g. from sharp-edged objects) are the main causes of corrosion.
- For sealing, choose appropriate seals.
- For flange mounting, only use sealings with a sufficiently large inner diameter and centre them; contact with the diaphragm leads to measuring deviations.
- When using soft or PTFE sealings, observe the instructions of the sealing manufacturer, particularly with regard to tightening torque and load cycles.
- For installation, in accordance with the fitting and flange standards the appropriate fastenings, such as screws and nuts, must be used.
- Model PG43SA-S with vent hole: The vent hole at the underside of the case must not be closed.
- Nominal position per EN 837-3 / 9.6.6 figure 7: 90° (⊥) (only applies to model PG43SA-S)

5.3 Mounting instructions for instruments with EHEDG and 3-A

Observe the following instructions, especially for EHEDG-certified and 3-A-compliant instruments.

- To maintain the EHEDG certification, one of the EHEDG-recommended process connections must be used. These are marked with the logo in the data sheet.

5. Commissioning, operation

- To maintain the conformity to the 3-A standard, a 3-A marked process connection must be used. These are marked with the logo in the data sheet.
- Mount the instrument with minimal dead space and able to be cleaned easily.
- The mounting position of the instrument, welding socket and instrumentation T-piece should be designed to be self-draining.
- The mounting position must not form a draining point or cause a basin to be formed.
- With the process connection via an instrumentation T-piece, the branch L of the T-piece must not be longer than the diameter D of the T-piece ($L \leq D$). If the instrument protrudes into the dead leg (reaching into the branch), this diameter “d” shall be subtracted from the inner diameter of the branch D (rule: $L \leq D - d$).
- The welding sockets shall comply with EHEDG guidelines doc. 9 and doc. 35. With TIG welding, the surface roughness must be $Ra \leq 0.8 \mu\text{m}$. The weld seams must be free from defects on the wetted and non-wetted contact sides and must be butt welded.
- For tank-mounted instruments the tank cleaning devices shall be positioned in such a way that the instrument can be assessed and is perfectly cleaned.
- Tank-mounted instruments shall be positioned in such a way to ensure that the spray from the cleaning device can reach all surfaces.

Special instructions for in-line versions

- In-line access units (e.g. VARINLINE®) and instrumentation T-pieces in horizontal pipes should be installed with a slight pipe inclination for self-draining.
- With in-line access units (e.g. VARINLINE®), use the original seals from the manufacturer and observe their installation instructions.

5. Commissioning, operation

EN

5.4 Permissible ambient and operating temperatures

When mounting the pressure gauge it must be ensured that, taking into consideration the influence of convection and heat radiation, no deviation above or below the permissible ambient and media temperatures can occur. The influence of temperature on the indication accuracy must be observed.

In order to avoid any additional heating, the instruments must not be exposed to direct solar irradiation while in operation.

5.5 Permissible vibration load at the installation site

The instruments should always be installed in locations free from vibration.

5.6 Commissioning

During the commissioning process pressure surges must be avoided at all costs.

5.7 External zero point setting (only model PG43SA-S with vent hole)

Setting range $\pm 15^\circ \triangleq \pm 5.5\%$ of the span

External zero point setting is not available for the model PG43SA-S with vent valve.

The zero point setting can be made with a slotted screwdriver after the removal of the sealing plug on the top of the case.

5.8 Autoclaving

The following is valid only for instruments with the option “Autoclavable version”, identifiable through the dial marking “AUTOCLAVABLE”.

- Model PG43SA-S with vent valve: To carry out the autoclaving process correctly, the vent valve must be open. See chapter 2.3 “Model PG43SA-S with vent valve”.
- Model PG43SA-C with vent valve: Ensure that, before the start of the autoclaving process, the vent valve has been closed. See chapter 2.4 “Model PG43SA-C with vent valve”.
- Observe the permissible conditions of autoclaving, see chapter 9.3 “Permissible temperature”
- During the autoclaving process, no weight-loading must be exerted via the case and the instrument must not lie on its window.
- Model PG43SA-S with vent hole: In case of a zero offset, carry out a zero point setting. See chapter 5.7 “External zero point setting”.
- Model PG43SA-C with vent valve: Allow the instrument to cool to ambient temperature. By briefly opening the vent valve, the pressure difference between the atmosphere and the sealed case will be eliminated. See chapter 2.4 “Model PG43SA-C with vent valve”.
- After the autoclaving cycles, there may be a slight bulging of the window. With an increasing number of autoclaving cycles, whitish residues can build up on the window and pointer. Hairline cracks may appear in the dial paint. It is recommended to subject the measuring instrument to a visual inspection for external damage after each autoclaving cycle.

In case of external damage, the instrument has to be taken out of service, if necessary, since the 3-A conformity may be impaired.

It is recommended to carry out a regular functional test after 30 autoclaving cycles. For this the instrument must be disconnected from the process to check with a pressure testing device. If this functional test is not passed, then the manufacturer must be contacted.

6. Faults

Following any failure or shutdown of the plant all instruments should be checked and, if necessary, replaced before recommissioning the plant.

7. Maintenance and cleaning

7.1 Maintenance

The instrument is maintenance-free. Checks should be carried out on a regular basis to ensure the measurement accuracy of the instrument. The checks or calibrations must be carried out by qualified skilled personnel with the appropriate equipment.

- The indicator should be checked once or twice every year. For this the instrument must be disconnected from the process to check with a pressure testing device.
- Repairs should only be carried out by the manufacturer or, following prior consultation with the manufacturer, by correspondingly qualified skilled personnel.
- With contaminated, viscous or crystallising media, it may be necessary to clean the diaphragm and/or instrument from time to time, see chapter 7.2 “Cleaning”.

7. Maintenance and cleaning

7.2 Cleaning

EN

NOTE!

Damage to property due to improper cleaning

Improper cleaning may lead to damage to the instrument.

- ▶ Do not use any aggressive cleaning agents or cleaning agents that cause corrosion or that are abrasive for the wetted parts. Only use cleaning agents suitable for the materials used in the instrument.
- ▶ Do not use any hard or pointed objects for cleaning.
- ▶ Do not use any abrasive cloths or sponges.

Cleaning in place (CIP) cleaning process

The instrument does not have to be dismantled for cleaning. It can be cleaned in place together with the pipeline or tank.

Avoid thermal shocks or fast changes in the temperature. The temperature difference between the cleaning agent and rinsing with clear water should be as low as possible. Negative example: cleaning with 80 °C and rinsing at 4 °C with clear water.

Cleaning after dismantling the instrument

1. Before cleaning, correctly disconnect the instrument from the pressure supply.
2. Clean the outside of the instrument with a moist cloth. If cleaning from outside ("wash down") is permitted, observe the permissible temperature and ingress protection. Make sure that the lowest ingress protection of the instrument is taken into account.
3. With contaminated, viscous or crystallising media, it may be necessary to clean the diaphragm from time to time. Only ever remove deposits from the diaphragm with a soft brush and suitable cleaning agent.

For sterilisation in an autoclave, observe the instructions in chapter 5.8 "Autoclaving".

8. Dismounting, return and disposal

WARNING!

Physical injuries and damage to property and the environment through residual media

Residual media in the dismantled instrument can result in a risk to persons, the environment and equipment.

- ▶ Wash or clean the dismantled instrument, in order to protect persons and the environment from exposure to residual media.

EN

8.1 Dismounting

Only disconnect the measuring instrument once the system has been depressurised.

8.2 Return

Strictly observe the following when shipping the instrument:

- All instruments delivered to WIKA must be free from any kind of hazardous substances (acids, bases, solutions, etc.) and must therefore be cleaned before being returned, see chapter 7.2 "Cleaning".
- When returning the instrument, use the original packaging or a suitable transport packaging.

INFORMATION!

With hazardous substances, include the material safety data sheet for the corresponding medium.

8.3 Disposal

Incorrect disposal can put the environment at risk.

Dispose of instrument components and packaging materials in an environmentally compatible way and in accordance with the country-specific waste disposal regulations.

9. Specifications

9. Specifications

EN

9.1 Pressure limitation

Model	Steady	Fluctuating	
		Scale range < 4 bar	Scale range ≥ 4 bar
PG43SA-S	Full scale value	Full scale value	2/3 x full scale value
PG43SA-C	3/4 x full scale value	2/3 x full scale value	

9.2 Overload safety

Model	Overload safety ¹⁾
PG43SA-S	2 x full scale value, max. 40 bar or max. pressure rating (PN) of the process connection
PG43SA-C	5 x full scale value, max. 40 bar or max. pressure rating (PN) of the process connection

1) see PN in the tables of data sheets PM 04.15, PM 04.16

9.3 Vacuum resistance (only model PG43SA-S)

All scale ranges are suitable for the vacuums typical of CIP or SIP cleaning. Positive scale ranges < 6 bar [< 100 psi] are optionally available in a vacuum-resistant version (to -1 bar).

9.4 Permissible temperature

Place, situation	Permissible temperature
Ambient	-20 ... +60 °C [-4 ... +160 °F]
Medium	≤ 150 °C [≤ 302 °F]
CIP and SIP	150 °C [302 °F] continuously for wetted parts
Autoclaving ¹⁾	≤ 134 °C [≤ 273 °F], ≤ 20 minutes
Storage, transport	-40 ... +70 °C [-40 ... +158 °F]

1) only for autoclavable version

9. Specifications

EN

9.5 Temperature effect

Model	Deviation from the reference temperature (+20 °C) at the measuring system
PG43SA-S	$\leq \pm 0.5 \text{ \%}/10 \text{ K}$ of full scale value
PG43SA-C	$\leq \pm 0.8 \text{ \%}/10 \text{ K}$ of full scale value

9.6 Ingress protection per IEC/EN 60529

Model PG43SA-S with vent hole: IP54

Model PG43SA-S with closed vent valve: IP68

Model PG43SA-S with open vent valve: IP54

Model PG43SA-C without vent valve: IP66

Model PG43SA-C with closed vent valve: IP68

Model PG43SA-C with open vent valve: IP54

For further specifications, see WIKA data sheets PM 04.15, PM 04.16 and the order documentation.

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VARINLINE® is a registered trademark of the company GEA Tuchenhausen.
BioControl® is a registered trademark of the company NEUMO.
WIKA Instruments Ltd